

HK GREATMADE TECH LTD

Keep pace with your need



Contact: Shirley Wang
 E-mail : info@greatmade.com , sales@greatmade.com
 Website : www.greatmade.com , www.greatmade.divtrade.com
 Tel : [86-0755-28144774](tel:86-0755-28144774), Fax: [86-0755-27185453](tel:86-0755-27185453)
 Hotmail : greatmade@hotmail.com
 Skype : Greatmade01

ECG cable and Electrode (心电图机电缆及导联线)

Burdick

(伯迪克)

OEM P/N	Greatmade No.	Description (Cable and connector)
EK-10	EC001A	Burdick EKG cable , IEC , DB15M> DIN3.0
	EC001AD	Burdick EKG cable , AHA , DB15M> DIN 3.0
	EC001B	Burdick EKG cable , AHA , DB15M> Banana 4.0
	EC001BI	Burdick EKG cable , IEC , DB15M> Banana 4.0




Fukuda

(福田)

OEM P/N	Greatmade No.	Description (Cable and connector)
FX-101	EC002A	Fukuda EKG cable , IEC , DB15M> Din 3.0
	EC002AD	Fukuda EKG cable , AHA , DB15M> DIN 3.0
	EC002B	Fukuda EKG cable , AHA , DB15M> Banana 4.0
	EC002BI	Fukuda EKG cable , IEC , DB15M> Banana 4.0



Suitable for **FUKUDA** FCP-7411, FX-7402, FX-7202, FX-2111
 FX-2155, FX-3010, FX-7102



GE Marquette

(GE-马葵)

OEM P/N	Greatmade No.	Description (Cable and connector)
	EC003A	MAC-500/MAC-1200 Multi-Link Cable , AHA
	EC003AI	Mac-500/ MAC-1200 Multi-Link Cable , IEC
412680-001	EC004	Multi-Link Leadwire set , Banana ending ,AHA
	EC004I	Multi-LiNK Leadwire set , Banana ending,IEC



GE Marquette One piece (GE-马葵 一体式)

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
	EC005A	<u>MAC-500/MAC-1200</u> One-piece EKG CABLE and leadwires <u>AHA</u> , <u>Banana</u>
	EC005AI	<u>MAC-500/MAC-1200</u> One-piece EKG CABLE and leadwires <u>IEC</u> , <u>Banana</u>



Kanz (铃谦)

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
PC-104	EC006A	Kanz EKG cable , <u>IEC</u> , DDK-16PIN > <u>DIN 3.0</u>
	EC006AD	Kanz EKG cable , <u>AHA</u> , DDK-16PIN > <u>DIN 3.0</u>
	EC006B	Kanz EKG cable , <u>AHA</u> , DDK-16PIN > <u>Banana 4.0</u>
	EC006BI	Kanz EKG cable , <u>IEC</u> , DDK-16PIN > <u>Banana4.0</u>



Suitable for Kanz : 103, 106

Cardioline delta 1 plus , delta 3 plus , delta 60 plus .

Kanz (铃谦)

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
PC-109	EC0015A	Kanz EKG cable , <u>IEC</u> , DDK-16PIN > <u>Banana 4.0</u>
	EC015B	Kanz EKG cable , <u>AHA</u> , DDK-16PIN > <u>Banana 4.0</u>



Suitable with Kanz : PC-109

Nihon Kohden : (日本光电)

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
901D	EC007A	Nihon Kohden 10-lead EKG cable <u>IEC</u> , NK DB15M > <u>DIN3.0</u>
9130K	EC007B	NK9130k EKG cable , <u>IEC</u> , 12ft Grey & DB15M >> <u>Din3.0</u>
	EC007C	Nihon Kohden 10-lead EKG cable <u>AHA</u> , NK DB15M > <u>banana 4.0</u>
	EC007D	Nihon Kohden 10-lead EKG cable <u>IEC</u> , NK DB15M > <u>banana 4.0</u>
	EC007E	Nihon Kohden ECGLAB Traadmill test cable , IEC, DB15M- <u>Snap</u>
	EC007F	Nihon Kohden ECGLAB Traadmill test cable , AHA, DB15M- <u>Snap</u>
	EC007G	Nihon Kohden 10-Lead EKG cable , AHA , NK DB15M- <u>Grabber</u> .
	EC007H	Nihon Kohden 10-Lead EKG cable , IEC, NK DB15M- <u>Grabber</u> .

Compatibilities Use with Nihon Kohden ECG-9620, ECG-9010, ECG-9020, ECG-9022, ECG-9110, ECG-9130, ECG-9132, DonqJiang ECG-11A, ECG-11B, ECG-11D, ECG-32A, (II), ECG-120. BIOCARE ECG-101, ECG-101G, ECG-300, ECG-300G, ECG-1200, ECG-9803



NEC

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
YY-1063	EC008A	NEC EKG cable , <u>IEC</u> , Round-15pin > <u>DIN3.0</u>
	EC008B	NEC EKG cable <u>AHA</u> , Round-15pin > <u>Banana4.0</u>



Compatible with NEC: K-3321, 3330, 3370

(惠普 / 飞利浦)**Philips (HP)**

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
M3703A/C	EC009A	HP EKG cable , <u>IEC</u> , DB15m> <u>Banana 4.0</u>
	EC009B	HP EKG cable <u>AHA</u> , DB15m> <u>Banana 4.0</u>



Compatible with Philips(HP):1772A,M3703C M2462A

**(席勒)****Schiller**

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
AT-1	EC010A	schiller EKG cable , <u>IEC</u> , <u>DB15M</u> > <u>Banana 4.0</u>
	EC010B	schiller EKG cable <u>AHA</u> , <u>DB15M</u> > <u>Banana 4.0</u>
	EC010C	Schiller EKG Cable , <u>IEC</u> , <u>DB15M</u> > <u>Grabber</u>
	EC010D	Schiller EKG Cable , AHA, <u>DB15M</u> > <u>Grabber</u>



Suitable for Schiller Cardiovit : AT-102, CS-200, Cardiovit AT-10 plus , AT-104 PC, AT-110

**Petas**

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
	EC016A	Petas EKG cable , IEC, Round 12pin > Banana 4.0
	EC016B	Petas EKG cable <u>AHA</u> , Round 12pin > Banana 4.0

**ShangHai Kohden****(上海光电)**

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
	EC011A	ShangHai Kohden EKG Cable <u>IEC</u> , DB15M> <u>DIN 3.0</u>
	EC011B	ShangHai Kohden EKG Cable <u>AHA</u> , DB15M> <u>Banana 4.0</u>
	EC011C	ShangHai Kohden EKG Cable , AHA, DB15M> <u>DIN 3.0</u>
	EC011B	ShangHai Kohden EKG Cable ,IEC, DB15M> <u>Banana 4.0</u>

Compatible with ShaiHai Kohden : ECG-6511, ECG-6151, 8110P/K, 6353,



MORTARA

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
	EC019A	Petas EKG cable , IEC, DB15M > Banana 4.0 , 4.7KQ resistance
	EC019B	Petas EKG cable AHA , DB15M > Banana 4.0 , 4.7KQ resistance



(成人肢体夹)

Limb Electrode

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
	EC012A	Limb Electrode (4pcs/set) Perfect for 3.0 DIN
	EC012B	Limb Electrode (4pcs/set) Perfect for 4.0 Banana
	EC012C	Limb Electrode with Screw (4pcs/set),Suitable for 3.0 Din and 4.0 Banana



EC012A Suitable for all 3.0 Din Plug of EKG Cable

EC012B Suitable for All 4.0 Banana Plug of EKG Cable .

EC012C For Multi-use , Suitable for all 3.0DIN , 4.0 Banana , Snap and Grabber

Pediatric Limb Electrode


(儿童肢体夹)

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
	EC012D	Limb Electrode with Screw (4pcs/set),Suitable for 3.0 Din and 4.0 Banana



Suction Electrode (吸球电极)

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
	EC013A	Suction Electrodes (6 pcs/ Set), perfect for 3.0 DIN
	EC013B	Suction Electrodes (6 pcs/ Set), perfect for 4.0 Banana
	EC013C	Suction Electrodes with Screw (6 pcs/ set),Suitable for 3.0 Din and 4.0 Banana, Snap



EC013A Suitable for all 3.0 Din Plug of EKG Cable

EC013B Suitable for All 4.0 Banana Plug of EKG Cable .

EC013C For Multi-use , Suitable for all 3.0DIN , 4.0 Banana , Snap and Grabber

EKG Cable Electrode Adapter

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
	EC017A	Electrode Adapter , <i>3.0DIN-> Snap , L=10CM</i>
	EC017B	Electrode Adapter <i>4.0 Banana -> Snap , L=10CM</i>
	EC017C	Electrode Adapter , <i>3.0DIN->Grabber , L=10CM</i>
	EC017D	Electrode Adapter <i>4.0 Banana-> Grabber , L=10CM</i>



Disposable Electrode (一次性电极片)

OEM P/N	Greatmade No.	Desceiption (Cable and connector)
	EC014	Disposable Electrodes , 1000 pcs/ box

