

| Compatible Brand | Item Number | CONNECTOR | | Compatible Models Number | Discription | |
|------------------|-------------|---|---|---|--------------------|---------------------------|
| BCI | SA0002A |  |  | (Biochem International): 3401,3304,3303,3302,3301, 3300 | Blue PU cable/2.4M | DB9M>DB9F >>BCI |
| | | | | | grey PU cable/2.4M | |
| BCI | SA0002B |  |  | (Biochem International): 6100,9100,Advisor,Autocorr, Mini Torr BCI,Mini Torr Plus | Blue PU cable/2.4M | 7P>DB9F >>BCI |
| | | | | | grey PU cable/2.4M | |
| Biocare | SA0073 |  |  | Biocare(China)邦健: BM9000 Jumbor 精博 | Blue PU cable/2.4M | 6P>DB9F >>BCI |
| | | | | | grey PU cable/2.4M | |
| Bionet | SA0040A |  |  | BM3Vet,BM31 | Blue PU cable/2.4M | 7P>DB9F >>Bionet |
| | | | | | grey PU cable/2.4M | |
| Bionet | SA0040B |  |  | BM3,BM3 Plus | Blue PU cable/2.4M | 圆6P>DB9F >>Bionet |
| | | | | | grey PU cable/2.4M | |
| Biosys | SA0039 |  |  | BPM300 | Blue PU cable/2.4M | MD6M>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |
| Bruker | SA0050 |  |  | Bruker all | Blue PU cable/2.4M | AMP14P>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |
| Colin | SA0029 |  |  | Datascope:Accutorr Plus, Nellcor:NPB-290/5,N-390,N-3000, N-6000, Nihon Kohden:BSM-410A, BSM-4114A, Colin BP88/306 | Blue PU cable/2.4M | AMP14P>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |
| Creative | SA0043 |  |  | Creative:PC9000,UP6000, TOB-MED 同科博越:TOM MS-6-P, EDAN 理邦: M9 BIOCARE:BM9000, Goldway CONTEC | Blue PU cable/2.4M | 6P>DB9F >>BCI |
| | | | | | grey PU cable/2.4M | |
| CSI | SA0011B |  |  | CSI(Criticare Systems International):506 Spot CSI, 507SD, 503 Spot | Blue PU cable/2.4M | DB9M>DB9F >>CSI |
| | | | | | grey PU cable/2.4M | |
| CSI | SA0011A |  |  | CSI(Criticare Systems International):Scholar,Poet(602- 3A),1100,581,508,507, 507SD, 504,503,501, Mindray 500A | Blue PU cable/2.4M | 1B5P>DB9F >>CSI |
| | | | | | grey PU cable/2.4M | |
| Datascope | SA0012 |  |  | Datascope:Accuat, Accutorr3/4 SAT Accutorr Plus Passport(Database Spo2) | Blue PU cable/2.4M | DIN8P>DB9F >>Datascope |
| | | | | | grey PU cable/2.4M | |
| Datex | SA0003 |  |  | Datex-Oheda:AS/3,CS/3,CH/S, CH/2S,CH-RS,CG/2S,Cardiopl, CardioplII,Cardiopl5, Capnomac Uitima,Light Monitor, M-Estpr,OscarII,OSP-200 Satlite Trans,Satlite Trans Plus, ULTSVI,ULTSO,ULTXIO, ULTS00 | Blue PU cable/2.4M | 10P>DB9F >>Datex |
| | | | | | grey PU cable/2.4M | |
| | SA0003-OXY |  |  | | Blue PU cable/2.4M | 10P>8J >>Datex |
| | | | | | grey PU cable/2.4M | |

| Compatible Brand | Item Number | CONNECTOR | | Compatible Models Number | Discription | |
|------------------|-------------|---|---|--|--------------------|----------------------------------|
| Datex -OHMEDA | SA0009 |  |  | Suitable for Datex-Ohmeda:Ohmeda 3700,Ohmeda3710,Ohmeda 3740, Ohmeda 3775,Biox 3740,Biox 3760 3800OxiCap,4700OxiCap,5250RGM, Rascal II GE Medical Systems:Corometrics 511, Mennen:Horizon XL (Ohmeda SpO2) | Blue PU cable/2.4M | Hyc>DB9F >>Ohmeda |
| | | | | | grey PU cable/2.4M | |
| | SA0009B |  |  | | Blue PU cable/2.4M | Hyc>8J >>Ohmeda |
| | | | | | grey PU cable/2.4M | |
| GE-Ohmeda | SA0009B-8J |  |  | Sensor solutions for GE Ohmeda TruSat pulse oximeter | Blue PU cable/2.4M | DB9M>8J >>Ohmeda |
| | | | | | grey PU cable/2.4M | |
| Dolphin(osi) | SA0036 |  |  | Dolphin ONE Oximeters | Blue PU cable/2.4M | AMP14P>DB9F >>Dolphin |
| | | | | | grey PU cable/2.4M | |
| EDNA | SA0075 |  |  | EDAN:M9 WANDONGÍ 4 4 BYS-8000£-MINDRAYÁ 4 4 PM-8000 | Blue PU cable/2.4M | 6P>DB9F >>BCI |
| | | | | | grey PU cable/2.4M | |
| GE Medical | SA0006 |  |  | GE Medical Systems:Eagle GE Medical Systems:Dash GE, Medical Systems:SOLAR GE, Medical Systems:TRAM Modular*50SL Series | Blue PU cable/2.4M | 11P>DB9F >>Nelcor |
| | | | | | grey PU cable/2.4M | |
| GE Medical | SA0006M |  |  | GE Medical Information technology:Oximax module | Blue PU cable/2.4M | 11P>DB9F >>Nelcor Oximax |
| | | | | | grey PU cable/2.4M | |
| GE Medical | SA0006-MB |  |  | GE Medical Systems:Dash, SOLAR Serie(Masimo module),İ İ -original or refurbished or repair only | Blue PU cable/2.4M | 11P>DB9F >>Masimo LNOP |
| | | | | | grey PU cable/2.4M | |
| | SA0006-MN |  |  | | Blue PU cable/2.4M | 11P>DB9F >>Nelcor |
| | | | | | grey PU cable/2.4M | |
| GE Ohmeda | SA0006-OXY |  |  | Sensor solutions for Eatex-Ohmeda S/5 FM ports: | Blue PU cable/2.4M | 11P>8J >>Ohmeda |
| | | | | | grey PU cable/2.4M | |
| Goldway | SA0015B |  |  | UT4000A,UT4000C,UT4000E, UT4000F, UT4000Arpo, UT4000Fpro | Blue PU cable/2.4M | 1B 5P>DB9F >>BCI |
| | | | | | grey PU cable/2.4M | |
| Goldway | SA0015 |  |  | Goldway:UT4000F, PACE TECH | Blue PU cable/2.4M | Redel 5P>DB9F >>BCI |
| | | | | | grey PU cable/2.4M | |
| Goldway | SA0015C |  |  | Goldway:monitor GOLDWAY MODEL UT6000A | Blue PU cable/2.4M | Redel 7P>DB9F >>Nelcor oximax |
| | | | | | grey PU cable/2.4M | |
| Invivo | SA0013 |  |  | 4500, 4500Plus, Millenia, Omega1445, Omni-TrakTVS, Sat Scout, original cable number 9383 | Blue PU cable/2.4M | 1B 7P>DB9F >>Nelcor |
| | | | | | grey PU cable/2.4M | |

| Compatible Brand | Item Number | CONNECTOR | | Compatible Models Number | Discription | |
|------------------|-------------|---|---|--|--------------------|-----------------------------------|
| Invivo | SA0013-MB |  |  | M8 original or refurbished or repair only | Blue PU cable/2.4M | 1B 7P>DB9F >>Masimo LNOP |
| | SA0013-MN |  |  | | Blue PU cable/2.4M | |
| Kontron | SA0031 |  |  | Kontron:7840 Kontron:7845 Kontron:Kolormon Kontron:Micromon Kontron:Minimon 7137 Kontron:Supermon Kontron:Trakmon | Blue PU cable/2.4M | 圆12P>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |
| Masimo | SA0016B |  |  | Masim all MDE:EscortII OPT30:Prism, GE Marquette:Pro1000, V300, INVIVO:T8 Schiller: ZOLL; Datascop: Trio, Passport2, Spectrum, original or refurbished or repair only | Blue PU cable/2.4M | AMP14P>6J >>Masimo LNOP |
| | | | | | grey PU cable/2.4M | |
| | SA0016A |  |  | | Blue PU cable/2.4M | AMP14P>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |
| MEK | SA0027 |  |  | MP-100, MP110, MP400, MP500, MP600, Mp1000 | Blue PU cable/2.4M | DB9M>DB9F >>MEK |
| | | | | | grey PU cable/2.4M | |
| MEK | SA0027B |  |  | MEK:NEW | Blue PU cable/2.4M | 圆8P>DB9F >>MEK |
| | | | | | grey PU cable/2.4M | |
| Mennen | SA0033B |  |  | | Blue PU cable/2.4M | 圆8P>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |
| Mennen | SA0033 |  |  | Horizon2000 | Blue PU cable/2.4M | 2B 8P>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |
| Mennen | SA0033 |  |  | 限维修, original or refurbished or repair only | Blue PU cable/2.4M | 13P>6J >>Masimo |
| | | | | | grey PU cable/2.4M | |
| Mindray | SA0018A |  |  | MEC-509,509+,503,507, OEM No.512A-30-06074 | Blue PU cable/2.4M | 1B 5P>DB9F >>CSI |
| | | | | | grey PU cable/2.4M | |
| Mindray | SA0018C |  |  | PM7000,6201,8000,9000, OEM NO.0010-21-11957 | Blue PU cable/2.4M | 1B 40 degree 6P>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |
| Mindray | SA0018B |  |  | PM5000,PM6000,PM8100 OEM No.0010-30-42737 | Blue PU cable/2.4M | 圆12P>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |
| Mindray | SA0018D |  |  | Be Compatible with Mindray Beneview T5/T6/T8. With Nellcor Technology,PM- 6800 OEM No.0010-30-42737 | Blue PU cable/2.4M | 圆7P>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |

| Compatible Brand | Item Number | CONNECTOR | | Compatible Models Number | Discription | |
|------------------|-------------|---|---|---|--|--|
| Mindray | SA0018E |  |  | VS800,PM80002,9201,(PM6201,7000, 8000,M1K0,M2K0 since 2006) OEM No.0010-30-42737 | Blue PU cable/2.4M grey PU cable/2.4M | 60degree 6P>DB9F >>Masimo |
| Mindray | SA0018C-MN |  |  | PM-9000,9100, 9200, (PM6201, 7000 8000,M1K0,M2K0, before 2006) OEM NO.0010-30-42737 | Blue PU cable/2.4M grey PU cable/2.4M | 40degree 6P>DB9F >>Masimo |
| Mindray | SA0018B-MN |  |  | PM-6000(masimo module, 12pin) OEM No.0010-30-42738 | Blue PU cable/2.4M grey PU cable/2.4M | 圆12P>DB9F >>Masimo |
| Mindray | SA0018D-MN |  |  | PM-6800(masimo module, 7pin) OEM No.0010-30-42738 | Blue PU cable/2.4M grey PU cable/2.4M | 圆7P>DB9F >>Masimo |
| Mindray | SA0018C-MB |  |  | PM9000 express, PM9000,9100,9200 (masimo module, 7pin) OEM No.0010-30-42738, 限维修, original or refurbished or repair only | Blue PU cable/2.4M grey PU cable/2.4M | 40 degree 6P>6J >>Masimo |
| Mindray | SA0018B-MB |  |  | PM-6000(masimo module, OEM No.0010-30-42738 限维修, original or refurbished or repair only | Blue PU cable/2.4M grey PU cable/2.4M | 圆12P>6J >>Masimo |
| Mindray | SA0018F |  |  | Use with PM-7000,PM-8000,PM-9000,MEC-1000,PM-9000Express,PM-600,new PM-6000, VS-800, for Nellcor oximax DS-100A spo2 sensor | Blue PU cable/2.4M grey PU cable/2.4M | 80degree 6P>DB9F >>NELLCOR OXIMAX |
| Nellcor | SA0001D |  |  | SCP-10. Datascope:Accutorr Pluss, Nellcor: NPB-290/5, N-390, N-395, N-3000,N-6000, NPB-290/5, N-390,N-395,N-3000,N-6000, Nihon Kohden:BSM-4104A, BSM-4114A,Colin BP88/306 | Blue PU cable/2.4M grey PU cable/2.4M | AMP14P>DB9F >>Nellcor |
| Nellcor | SA0001A |  |  | EC-8, Nellcor:N100,N200,N180,N-20,NPB- 40,NPB-75 | 蓝色PU/2.4M 灰色PU/2.4M | DB9M>DB9F >>Nellcor |
| Nellcor | SA0001E |  | OXIMAX adapter | Nellcor: N550 N-590 GE: Dinamap PRO(w/Oxi), Shiller All Nellcor & Nellcor equipped OEM nonitors not listed above | Blue PU cable/2.4M grey PU cable/2.4M | DB9M 7Pin>DB9F >>Nellcor old Nellcor to Oxima Tech Oximeter |
| Nellcor | SA0001B |  |  | DEC-8, All Nellcor & Nellcor equipped OEM nonitors not listed above | Blue PU cable/2.4M grey PU cable/2.4M | DB9M 9Pin>DB9F >>Nellcor |
| Nellcor | SA0001C |  |  | Nellcor: N550,N-595, GE: Dinamap PRD(W/Oxi),Schiller | Blue PU cable/2.4M grey PU cable/2.4M | AMP14P>DB9F >>Nellcor oxinax |
| Nellcor | SA0001F |  |  | Nellcor Pre-amp spo2 adapter cable , Lemo 12pin to DB9 male for N100 (w High Sat),N200 | Blue PU cable/2.4M grey PU cable/2.4M | Lemo12pin>DB9F |
| Neusoft | SA0066 |  |  | NSC series patient monitor | Blue PU cable/2.4M grey PU cable/2.4M | 圆12P>DB9F >>Nellcor |

| Compatible Brand | Item Number | CONNECTOR | | Compatible Models Number | Discription | |
|------------------|-------------|---|--|---|--------------------|------------------------------|
| Nihon Konden | SA0019 |  |  | AL-800PL Module,BSM-1101; LifeScope EC, BSM-2101/2 LifeScope L, BSM-3101 LifeScope LC,BSM-8200/1 LifeScope12, BSM83-01 /2 LifeScope9, SM-8800 LifeScope14,OLV-1100A | Blue PU cable/2.4M | 圆10P>DB9F >>Nihon Konden |
| | | | | | grey PU cable/2.4M | |
| Nihon Konden | SA0019C |  |  | Bedside Monitor:BSM4101,4103,4111, 4113 Life scope P, BSM-5105,5135 Life scope A, BSM-2351,2353 Life scope L,BSM-2301,2303 Life scope I | Blue PU cable/2.4M | 圆10P>DB9F >>Nihon Konden |
| | | | | | grey PU cable/2.4M | |
| Nihon Konden | SA0019B |  |  | Nihon Kohden JL-302T spo2 adapter cable | Blue PU cable/2.4M | 方14P>DB9F >>Nihon Konden |
| | | | | | grey PU cable/2.4M | |
| Nihon Konden | SA0019D |  |  | 8604,8604D | Blue PU cable/2.4M | 3M 20pin >Db9 |
| | | | | | grey PU cable/2.4M | |
| Nonin | SA0020A |  |  | 8600 series, 8500,8700,8800,9700 | Blue PU cable/2.4M | 1B 6P>DB9F >>Nonin |
| | | | | | grey PU cable/2.4M | |
| Novamatrix | SA0008 |  |  | Novamatrix:505,510,511,515B/C 520,7000,7100,Oxypleth; Spacelabs Medical:90351-D/G,90465,90466,90467,90489,90651-07/08,90419-01,90651A-07 | Blue PU cable/2.4M | HYP7P>DB9F >>Novamatrix |
| | | | | | grey PU cable/2.4M | |
| Novamatrix | SA0008B |  |  Hyp母座 | Oxypleth(Extension Cable),7100(Extension Cable), 7000(Extension Cable),520(Extension Cable),515B/C(Extension Cable),505(Extention Cable),500 | Blue PU cable/2.4M | HYP7P>HYP7J >>Novamatrix |
| | | | | | grey PU cable/2.4M | |
| Pace Tech | SA0030C |  |  | Pace Tech:All w/black face(except Vitalmax 800 w/Nellcor SpO2), Generra:710,200,300,400,810,530, | Blue PU cable/2.4M | 5P>DB9F >>BCI |
| | | | | | grey PU cable/2.4M | |
| Pace Tech | SA0030B |  |  | | Blue PU cable/2.4M | 7P>DB9F >>BCI |
| | | | | | grey PU cable/2.4M | |
| Philips | SA0005C |  |  | Philips Medical:M2601A,M3000A,M3500 B,M4735A | Blue PU cable/2.4M | 8PIN >DB9 Female>>Nellcor |
| | | | | | grey PU cable/2.4M | |
| Philips | SA0005C-M |  |  | With Masimo Technology Use with Masimo Spo2 sensor LNCS DC-I Af074 | Blue PU cable/2.4M | 8PIN >DB9 Female>>masimo |
| | | | | | grey PU cable/2.4M | |
| Philips | SA0005D |  |  | Philips Medical:M2601A,M3000A,M3500 B,M4735A | Blue PU cable/2.4M | 8P>8J >>Philips |
| | | | | | grey PU cable/2.4M | |
| Philips | SA0005B |  |  | Philips Medical:78352C,78354C,78834C, M1020A,M1025B,M1350B,M1722 A/B, M1732A/B,M2475B, Viridia A1, Viridia A3 | Blue PU cable/2.4M | 圆12P>DB9F >>Nellcor |
| | | | | | grey PU cable/2.4M | |