



SURGICAL STAPLER

Titan E™

Disposable Endoscopic Powered Cutter Stapler and Cartridge

Applied to tissue transection, resection and anastomosis in thoracic surgery, abdominal surgery and gynecologic endoscopy.



Reloads

Model	Color	Effective Length	Staple Hight	Staple Qty	Trocar	Model	Length
ZEQH 45W	White	48.5	2.5	70	12	ZEQHB 45W	253
ZEQH 45B	Blue		3.5			ZEQHL 45W	353
ZEQH 45D	Gold		3.8			ZEQHB 60W	253
ZEQH 45G	Green		4.1			ZEQHL 60W	353
ZEQH 45W	White	60	2.5	88			
ZEQH 45B	Blue		3.5				
ZEQH 45D	Gold		3.8				
ZEQH 45G	Green		4.1				

Stapler

Vega E™

Disposable Endoscopic Powered Cutter Stapler and Cartridge

Used in the tissue transection, resection and anastomosis of endoscopic surgery, gynaecology, paediatric and chest surgery.



Reloads

Model	Color	Suture Length	Cutting Length	Staple Height	Trocar	Model	Type	Length
ZEQN 30KX	White	30	30	2.0/2.5/3.0	12	ZEQN SW	Short	277
ZEQN 30JX	Purple			3.0/3.5/4.0				
ZEQN 45KX	White	45	45	2.0/2.5/3.0	15	ZEQN MW	Standard	377
ZEQN 45JX	Purple			3.0/3.5/4.0				
ZEQN 45GX	Black	60	60	4.0/4.5/5.0	15			
ZEQN 60KX	White			2.0/2.5/3.0	12	ZEQN LW	Extra-long	477
ZEQN 60JX	Purple			3.0/3.5/4.0				
ZEQN 60GX	Black	45	45	4.0/4.5/5.0	15			
ZEQN 45AX	White			2.5	12			
ZEQN 45BX	Blue			3.5				
ZEQN 45CX	Gold			4.0	15			
ZEQN 45DX	Green	60	60	4.8				
ZEQN 60AX	White			2.5	12			
ZEQN 60BX	Blue			3.5				
ZEQN 60CX	Gold			4.0	15			
ZEQN 60DX	Green			4.8				

Titan™

Disposable Endoscopic Cutter Stapler and Cartridge

Applied to tissue transection, resection and anastomosis in thoracic surgery, abdominal surgery and gynecologic endoscopy.



Reloads

Model	Color	Effective Length	Staple Hight	Staple Qty	Trocar	Model	Length
ZQH 45W	White	48.5	2.5	70	12	ZQHB 45W	516
ZQH 45B	Blue		3.5			ZQHL 45W	616
ZQH 45D	Gold		3.8			ZQHB 60W	528
ZQH 45G	Green		4.1			ZQHL 60W	628
ZQH 45W	White	60	2.5	88			
ZQH 45B	Blue		3.5				
ZQH 45D	Gold		3.8				
ZQH 45G	Green		4.1				

Stapler

Vega™

Disposable Endoscopic Cutter Stapler and Cartridge

Used in the tissue transection, resection and anastomosis of endoscopic surgery, gynaecology, paediatric and chest surgery.



Reloads

Model	Color	Suture Length	Cutting Length	Staple Hight	Trocar	Model	Type	Length
ZQN 30KX	White	30	30	2.0/2.5/3.0	12	ZQN SW	Short	277
ZQN 30JX	Purple			3.0/3.5/4.0				
ZQN 45KX	White	45	45	2.0/2.5/3.0		ZQN MW	Standard	377
ZQN 45JX	Purple			3.0/3.5/4.0				
ZQN 45GX	Black			4.0/4.5/5.0	15			
ZQN 60KX	White	60	60	2.0/2.5/3.0	12	ZQN LW	Extra-long	477
ZQN 60JX	Purple			3.0/3.5/4.0	15			
ZQN 60GX	Black			4.0/4.5/5.0				
ZQN 45AX	White	45	45	2.5	12			
ZQN 45BX	Blue			3.5				
ZQN 45CX	Gold			4.0				
ZQN 45DX	Green			4.8	15			
ZQN 60AX	White	60	60	2.5	12			
ZQN 60BX	Blue			3.5				
ZQN 60CX	Gold			4.0				
ZQN 60DX	Green			4.8	15			

Disposable Circular Stapler

Applied in digestive tract reconstructive surgery for end-to-end, end-to-side and side-to-side anastomosis.



Model	Outer Diameter	Cutting Diameter	Staple Height	Staple Width	Staple Qty
ZGT 21C	21.2	12.4	4.5	3.7	27
ZGT 24C	24.1	15.0	4.5	3.4	33
ZGT 25C	25.1	16.4	4.5	3.8	33
ZGT 26C	26.1	17.0	4.5	3.8	33
ZGT 29D	29.6	20.3	4.8	3.9	36
ZGT 31D	31.0	21.0	4.8	3.9	36
ZGT 32E	32.1	22.0	5.0	3.9	42

Disposable Circular Stapler

Applicated in the whole digestive tract for end-to-end, end-to-side anastomose.



Model	Outer Diameter	Cutting Diameter	Staple Height	Staple Width	Staple Qty
ZGX 17C	17.6	10.0	4.5	3.2	16
ZGX 19C	19.6	11.8	4.5	4.0	16
ZGX 21C	21.6	12.4	4.5	4.0	16
ZGX 24C	24.3	15.0	4.5	4.0	18
ZGX 25C	25.5	16.4	4.5	4.0	20
ZGX 26C	26.0	17.0	4.5	4.0	20
ZGX 29D	29.7	20.7	4.8	4.0	24
ZGX 32E	31.9	23.1	5.0	4.0	26
ZGX 34E	34.0	24.0	5.0	4.0	32

Disposable PPH Stapler

Applied to anal hemorrhoidectomy and excision of rectal prolapse.



Model	Outer Diameter	Cutting Diameter	Staple Height	Staple Width	Staple Qty
ZPPH III 34	33.5	23.1	3.8	3.8	48
ZPPH III 36	35.0	24.8	4.4	3.9	48

Disposable PPH Stapler

Suitable for anal hemorrhoidectomy, rectal prolapse,etc



Model	Outer Diameter	Cutting Diameter	Staple Height	Staple Width	Staple Qty
ZPPH II 32	32.5	24.1	3.8	4	32
ZPPH II 34	33.5	25.6	4.2	4	32

Disposable Linear Cutter Stapler and Cartridge

Applied to dissection, resection and suture of the digestive tract tissues of the stomach and intestines.



Reloads				Stapler			
Model	Color	Staple Height	Staple Qty	Model	Basic Size	Cutting Length	Closure Height
ZZQSR 60C	Blue	3.8	64	ZZQS 60C	255	60	1.5
ZZQSR 60D	Green	4.8	4	ZZQS 60D			2.0
ZZQSR 80C	Blue	3.8	84	ZZQS 80C	278	80	1.5
ZZQSR 80D	Green	4.8	84	ZZQS 80D			2.0
ZZQSR 100C	Blue	3.8	104	ZZQS 100C	315	100	1.5
ZZQSR 100D	Green	4.8	104	ZZQS 100D			2.0

Disposable Auto Linear Stapler

Applied to the closure of the stump or incision of digestive tract.



Reloads				Stapler	
Model	Effective Length	Staple Height	Staple Qty	Model	
ZAZXR 30	31.5	4.5	11	ZAZX 30	
ZAZXR 45	42.5		15	ZAZX 45	
ZAZXR 60	59.0		21	ZAZX 60	
ZAZXR 90	86.5		31	ZAZX 90	

Disposable Curved Stapler and Cartridge (30)

Applied to anastomosis and resection of anal sphincter organizations and mucosa part.



Reloads

Model	Effective Length	Staple Height	Staple Weight	Staple Qty	Model
ZHXR 30	20	4.6	4	21	ZHX 30

Stapler

Disposable Curved Stapler and Cartridge (45)

Suitable for low anterior resection of rectal surgery.



Reloads

Model	Effective Length	Staple Height	Staple Weight	Staple Qty	Model
ZHXR 45	40	4.6	4	46	ZHX 45

Stapler

Disposable Trocar

Suitable for puncture of laparoscopic surgery, establishment of equipment access.



Model	External Diameter	Trocar Diameter	Trocar Length	Trocar Length	Total Length
ZTB5-100	5.5	6	100	140	170
ZTB10-100	10.3	11			180
ZTB12-100	12.8	13			
ZTB15-100	15.3	16			

Disposable Skin Stapler

Applied to the closure of the stump or incision of digestive tract.



Model	Staple Qty
ZPF 35W	35



Cooperative · Innovative · Excellent

ZENER MEDTEC (CHANGZHOU) CO., LTD
Bldg 9D, Jintong International Industrial Park, No. 8 Xihu
Road, Wujin Hi-Tech Industry Zone, Changzhou, China
info@zenermed.com 86-0519-83800520
www.zenermed.com