



Medplus-Life plus,your best choice

**Catalogue of NIBP Connector
Shenzhen Medplus Accessory Co.,Limited**





Contact:Tina Wu,E-mail:info@spo2sale.com.Tel:+86-755-29093174

Address:2nd Floor,2nd Building, 8th Industrial Park,
TianLiao Villiage, GongMing Town, Shenzhen Of China





NO.	P/N	Picture	OEM P/N	Description	
1	MPJ01			Locking male connector , Metal	
2	MPJ02		30068 (Connect with MPJ04); WA:5082-168	Female Luer Connector (For Nihon Kohden ,CSI , Datascope Air Hose)	
3	MPJ03		Connect with MPJ12	Male Slip Metal connector	
4	MPJ04		30065 (Connect with MPJ02); WA:5082-165	BCI Monitor (Old) Datascope(Old)	
5	MPJ05		300619 (Connect with MPJ13) WA:5082-161	Female screw connector ,Metal	
6	MPJ06				
7	MPJ07		330060(Connect with MPJ08)	Schiller Air hose connector	
8	MPJ08		33059 (Connect with MPJ07) WA:5082-179	Criticare, Datascope Drage r /Siemens Physio control Mindray, Philips/HP Spacelabs, Nellcor	
9	MPJ09			Female Bayonet, L=38mm	
10	MPJ10		Connect with MPJ11	Colin male luer Connector	
11	MPJ11		Connect with MPJ10	Colin female luer Connector	

NO.	P/N	Picture	OEM P/N	Description	
12	MPJ12		300669 (Connect with MPJ03) WA:5082-169	Female slip Metal connector	
13	MPJ13		300664 (Connect with MPJ05) WA:5082-164	Male screw connector, Plastic	
14	MPJ14		Connect with MPJ07; WA:5082-184	Bayonet male Plastic	
15	MPJ15		Connect with MPJ12 WA:5082-176		
15	MPJ16				
17	MPJ17		330091 (Connect with MPJ32); WA:5082-178	Nihon Kohden Datex(1), GE Solar, Dash, Eagle, TRAM	
18	MPJ18		330057	Plastic Female locking socket	
19	MPJ19			Redel male connector for Creative NIBP air hose	
20	MPJ20			Redel female connector for NIBP socket	
21	MPJ21		Connect with MPJ33	GE;AS/3, CS/3, S/5 Cardiocap 5/II, Light Eagle4000, Solar Tram1/2/300	
22	MPJ22			GE Solar, Dash, Eagle Series	
23	MPJ23			GE-Dash, Eagle, SOLAR Tram400/450/451 Tram800/850/851, Coro 120	
24	MPJ24		330073	GE-Dinamap, Pro, Compact, MPS Pro Care series	

NO.	P/N	Picture	OEM P/N	Description	
25	MPJ25			Philips/HP	
26	MPJ26		300829	Y adaptor-Adapter, Y, 2 TB to 1 TB, Unfi tted	
27	MPJ27			Drager/Siemens	
28	MPJ28			Two end thread	
29	MPJ29			GE Critikon Dinamap 8100 air hose connector	
30	MPJ30			Double tube hose internal connector	
31	MPJ31			Nihon Kohden(New)	
32	MPJ32		330064 (Connect with MPJ17) WA:5082-172	Zoll Series M Datex(1) Marquette(2)	
33	MPJ33		Connect with MPJ21	Pace Tech	
34	MPJ34		300668 (Connect with MPJ35) WA:5082-168	Mennen (Old) Schiller Bruker/Odam Spacelabs(Old)	
35	MPJ35		300665 (Connect with MPJ34) WA:5082-165	BCI Monitor (Old) Datascope(Old)	
36	MPJ36		Connect with MPJ35	Disposable Cuff	
37	MPJ37		330046	Adapter, Y, 2 Female Screw to 1 Male Bayonet	

NO.	P/N	Picture	OEM P/N	Description	
38	MPJ38		330083	Quick Connect Adapter, Y, 2 Male Submin to 1 Male Submin	
39	MPJ39		330082	Quick Connect Adapter, Y, 2 Male Submin to 1 Male Bayonet to 1 Female Submin	
40	MPJ40		330080	Quick Connect Adapter, Y, 2 Male Submin to 1 Male Bayonet	
41	MPJ41		330066	Adapter, Y, 2 Male Slip Luer to 1 Female Submin	
42	MPJ42		330049	Adapter, Y, 2 Female Slip Luer to 1 Male Slip Luer	
43	MPJ43		330069	Adapter, Y, 2 Male Slip Luer to 1 Female Submin 1 Male Submin	
44	MPJ44		330048	Adapter, Y, 2 Female Screw to 1 Unfi tted	
45	MPJ45		330084	Quick Connect Adapter, Y, 2 Male Submin to 1 Unfi tted	
46	MPJ46		330085	Quick Connect Adapter, X, 2 Male Submin to Submin Mated Pair	
47	MPJ47		330086	Quick Connect Adapter, X, 2 Male Submin to 2 Male Bayonet	
48	MPJ48		330087	Quick Connect Adapter, X, 2 Male Submin to 2 Screw Connectors	
49	MPJ49		300939	Adapter, X, 2 Female Submin to 2 Male Slip Luer	
50	MPJ50		330063	Adapter, X, 2 Male Slip Luer to Submin Mated Pair	

NO.	P/N	Picture	OEM P/N	Description	
51	MPJ51		330067	Adapter, X, 2 Male Slip Luer to 2 Male Bayonet	
52	MPJ52		300959	Tubing, Gray 2 in., 5/32 in. ID, Unfi tted	
53	MPJ53			Two-Tubing, Gray, Inner diameter: 5/32 in. Out diameter: 6mm*2, Purchasing Unit: M- Material: TPU	
54	MPJ54		300669	5/32 in. ID Barbed to Female Slip Luer, Metal	
55	MPJ55		300667	5/32 in. ID Barbed to Male Slip Luer, Metal,	
56	MPJ56		330092	Quick connect submin female 5/32"	
57	MPJ57			sphygmometer ball with valve	
58	MPJ58			Pressure gauge Pocket 50MM Scale Of sphygmometer	
59	MPJ59			sphygmomanometer bulb	
58	MPJ60			Dixtal male NIBP Connector	
59	MPJ61			Dixtal Female NIBP Connector	
58	MPJ62		Connect with MPJ63	Drager/Siemens	
59	MPJ63		Connect with MPJ62	Drager/Siemens	

NO.	P/N	Picture	OEM P/N	Description	
58	MPJ64			Female Bayonet,L=29mm	
59	MPJ65		Connect with MPJ66)		
58	MPJ66		Connect with MPJ650		
58	MPJ67		300664(Connect with MPJ05) WA:5082-164 L=22.6mm	Male screw Connector,metal	