



Electrode for Pelvic Floor Muscle Simulator

Description: The Caredmed EPFM electrodes are designed for pelvic floor muscle stimulation and biofeedback. It is supplied with an indicator stick which gives feedback for pelvic floor contractions whilst using it with a stimulation or EMG device and also allows it to be used as a Pelvic Floor Educator which means it can then also be used as a pelvic floor educator when practicing muscle contractions when not connected to the stimulation or feedback device. Single patient use only. This electrode's structure kindly meets human engineering and all materials are compliant with biocompatibility, it can be compatible with major brand of the market. Base on the intend application, the electrodes are classified as vaginal electrode, rectal electrode, vaginal pressure electrode and so on.



Product Part Number List

Item	P/N	Photo	Description
1	EPFM-V1		Vaginal electrode, smaller contact area
2	EPFM-V2		Vaginal electrode, larger contact area

3	EPFM-PV1		Vaginal pressure electrode
4	EPFM-R1		Rectal electrode, smaller contact
5	EPFM-R2		Rectal electrode with adjustable collar.

Cleaning

Caremed EPFM is designed for single person use and should not be autoclaved. It is manufactured from medical grade materials but not in guaranteed sterile conditions and we therefore recommend washing with soap and water prior to use.

Package

Initial pack: 1pcs/PE bag; outer package: 1pcs/white box;

Warranty

1 year

Shenzhen Caremed Medical Technology Co., Ltd

Add: East side, 3/F, C building, Kelunte low-carbon industrial park, Gaofeng community, Dalang office, Shenzhen city, China

Tel:+86755-36560977 Fax: +86755-28186486 Email:cm003@szcaremed.com

www.szcaremed.com