





UpnMed®

Design with heart, use with confidence

# ECG / EKG Cable with Leadwire

## Sweet Heart

This time we redefined the ECG/EKG Cable with Leadwires

Adhering to the concept of practicality and innovation, we abandon the old and unchanged traditional style. A bold heart-shaped appearance is used as a design element, giving a more reassuring quality and a higher-end texture.

Latest 107 days 36 modification

27 conferences 13 procedures

Craftsmanship is rigorous, only to guard carefully

CE ISO 13485





**NEW**

Sweet Heart ECG Junction Box

**In-mold injection molding**  
Durable and Tensile


Compact and lightweigh looking, no drop pressure, no winding wire, more convenient for doctors and nurses.

**Appearance  
P K**

Old-fashioned and heavy looking, a severe pulling feeling arised during use, inconvnience to use.



The new mark CF and defibrillation CF are added, which are easy to identify and convience for Doctor and Nurse to choose the correct one.

**Mark  
P K**


There is no mark to identy defibrillation or non-defibrillation, which is not conducive to intuitive selection and easy to make mistakes during use.



Integral molding process, IPX6, it will not effect the use if the surface accidentally sprayed with water.

**Waterproof  
P K**

The assembly process has poor waterproof effect, it cannot be used or damage the machine if the surface is accidentally sprayed with water.



Small body with large capacity. Upgraded Internal circuit design, can withstand high voltage AC 1500V, leakage current in line with AAMI EC53 standard.

**Quality  
P K**

The internal circuit can withstand high voltage AC 1000V, which cannot meet the regulatory standards.