

























































































Catalogue of IBP cable


















Brand	Part No	Description	Picture		Compatibility	
Biolight	MPI021A	Utah connectotr			Biolight A8;V6 Vital Sign View A2/A3/A5/A6	
	MPI021B	BD connector				
	MPI021C	Abbott/Medix connector				
	MPI021D	Philips connector				
	MPI021E	Edward/Baxter connector				
	MPI021F	Medex connector				
	MPI021G	PVB connector				
	MPI021H	Medex Logical/MX960				
Chenwei	MPI009A	Utah connectotr			Chenwei	
	MPI009B	BD connector				
	MPI009C	Abbott/Medix connector				
	MPI009D	Philips connector				
	MPI009E	Edward/Baxter connector				
	MPI009F	Medex connector				
	MPI009G	PVB connector				
	MPI009H	Medex Logical/MX960				



















Brand	Part No	Description	Picture		Compatibility
Choice	MPI013A	Utah connectr			
	MPI013B	BD connector			
	MPI013C	Abbott/Medix connector			
	MPI013D	Philips connector			
	MPI013E	Edward/Baxter connector			
	MPI013F	Medex connector			
	MPI013G	PVB connector			
	MPI013H	Medex Logical/MX960			
Colin	MPI030A	Utah connectr			
	MPI030B	BD connector			
	MPI030C	Abbott/Medix connector			
	MPI030D	Philips connector			
	MPI030E	Edward/Baxter connector			
	MPI030F	Medex connector			
	MPI030G	PVB connector			
	MPI030H	Medex Logical/MX960			



















Brand	Part No	Description	Picture		Compatibility
Corometrics	MPI020A	Utah connectotr			Corometrics
	MPI020B	BD connector			
	MPI020C	Abbott/Medix connector			
	MPI020D	Philips connector			
	MPI020E	Edward/Baxter connector			
	MPI020F	Medex connector			
	MPI020G	PVB connector			
	MPI020H	Medex Logical/MX960			
Creative	MPI024A	Utah connectotr			Creative
	MPI024B	BD connector			
	MPI024C	Abbott/Medix connector			
	MPI024D	Philips connector			
	MPI024E	Edward/Baxter connector			
	MPI024F	Medex connector			
	MPI024G	PVB connector			
	MPI024H	Medex Logical/MX960			



















Brand	Part No	Description	Picture		Compatibility	
Criticare	MPI023A	Utah connectotr			Criticare:Model 2200;CSI Monitor Scholar 8100 OEM Part No.: 650- 299-58	
	MPI023B	BD connector				
	MPI023C	Abbott/Medix connector				
	MPI023D	Philips connector				
	MPI023E	Edward/Baxter connector				
	MPI023F	Medex connector				
	MPI023G	PVB connector				
	MPI023H	Medex Logical/MX960				
Datascope	MPI028A	Utah connectotr			Datascope Passport V, V12, V21	
	MPI028B	BD connector				
	MPI028C	Abbott/Medix connector				
	MPI028D	Philips connector				
	MPI028E	Edward/Baxter connector				
	MPI028F	Medex connector				
	MPI028G	PVB connector				
	MPI028H	Medex Logical/MX960				



















Brand	Part No	Description	Picture	Compatibility
Datascope	MPI016A	Utah connectr		Datascope: System 80/82/83/84 90/90T/95/98/97IABP 850/870/P2/P3/2000 2000A/2001/2001A 2002/2002A/2200I, 3000/Passport; American Optical 1131 Module (Alternate Conn.). Colin 408BP, 400 Series. CWE Research Lab Monitor, Hewlett Packard Sandborn INEX. GE Marquette 3100series, SMEC 1300M (Balloon Pump)
	MPI016B	BD connector		
	MPI016C	Abbott/Medix connector		
	MPI016D	Philips connector		
	MPI016E	Edward/Baxter connector		
	MPI016F	Medex connector		
	MPI016G	PVB connector		
	MPI016H	Medex Logical/MX960		
Datex	PI001A	Utah connectr		Airshield, Burdick, Kone; Air-Shields AS-100 (141 Module), AS-441. Burdick M500 Series M544 Module (M560/M565/M575 Mainframe)." Datex-Ohmeda" AS/3, AS/5, ESTP Module, CC-104, Cardio-Cap CM- 104/Cap I and II/Cap V." Diascope Traveler, Kone" M500 Series M544 Module (M560/M565/M575 Mainframe)." Puritan - Bennett PB 240/245/245P. Simon son & Weil AS100 Series, Athena System, 8000, Quadrascope.
	MPI001B	BD connector		
	MPI001C	Abbott/Medix connector		
	MPI001D	Philips connector		
	MPI001E	Edward/Baxter connector		
	MPI001F	Medex connector		
	MPI001G	PVB connector		
	MPI001H	Medex Logical/MX960		



















Brand	Part No	Description	Picture	Compatibility		
Drager	MPI012A	Utah connectotr			Drager PM8014 & PM8060;Cato, Cicero EM, Vitara Parameterbox PB 8800; PM8010; PM 8040,Sulla, UM3, UM3.1	
	MPI012B	BD connector				
	MPI012C	Abbott/Medix connector				
	MPI012D	Philips connector				
	MPI012E	Edward/Baxter connector				
	MPI012F	Medex connector				
	MPI012G	PVB connector				
	MPI012H	Medex Logical/MX960				
Fukuda	MPI010A	Utah connectotr		Fukuda Denshi:DYNASCOPE /Expert DS5100E/DS5300W DS5400		
	MPI010B	BD connector				
	MPI010C	Abbott/Medix connector				
	MPI010D	Philips connector				
	MPI010E	Edward/Baxter connector				
	MPI010F	Medex connector				
	MPI010G	PVB connector				
	MPI010H	Medex Logical/MX960				



















Brand	Part No	Description	Picture	Compatibility		
GE	MPI022A	Utah connectr			GE-7000/Bard	
	MPI022B	BD connector				
	MPI022C	Abbott/Medix connector				
	MPI022D	Philips connector				
	MPI022E	Edward/Baxter connector				
	MPI022F	Medex connector				
	MPI022G	PVB connector				
	MPI022H	Medex Logical/MX960				
Guoteng	MPI025A	Utah connectr			Guoteng	
	MPI025B	BD connector				
	MPI025C	Abbott/Medix connector				
	MPI025D	Philips connector				
	MPI025E	Edward/Baxter connector				
	MPI025F	Medex connector				
	MPI025G	PVB connector				
	MPI025H	Medex Logical/MX960				



















Brand	Part No	Description	Picture		Compatibility
Hellige	MPI029A	Utah connectotr			Hellige 5th generation monitor SMU 610/611/612/613/639, Hon eywell/PPG SMU&VICOM series,Vitatek
	MPI029B	BD connector			
	MPI029C	Abbott/Medix connector			
	MPI029D	Philips connector			
	MPI029E	Edward/Baxter connector			
	MPI029F	Medex connector			
	MPI029G	PVB connector			
	MPI029H	Medex Logical/MX960			
Kontron	MPI019A	Utah connectotr			Kontron:Minimon 7131-7136, Fetalmon 2A 2BFetalogic 4000 Fetasonde 5/5A/5B, K2000 Ballon Pump, KAAT II/Module 10-7286,Supermon
	MPI019B	BD connector			
	MPI019C	Abbott/Medix connector			
	MPI019D	Philips connector			
	MPI019E	Edward/Baxter connector			
	MPI019F	Medex connector			
	MPI019G	PVB connector			
	MPI019H	Medex Logical/MX960			



















Brand	Part No	Description	Picture		Compatibility	
Lohmeier	MPI007A	Utah connectr			M010	
	MPI007B	BD connector				
	MPI007C	Abbott/Medix connector				
	MPI007D	Philips connector				
	MPI007E	Edward/Baxter connector				
	MPI007F	Medex connector				
	MPI007G	PVB connector				
	MPI007H	Medex Logical/MX960				
Marquette	MPI002A	Utah connector (700078-001)			GE Solar,Dash,Eagle,Tram; Eagle/dash1000/2000/ 3000/3100/4000,Solar 7000/8000/9000,Tram Module 100-600,Cath-Lab, MAC-Lab Dash †/Tram(all versions),E-PSMP/E-PRESTN/PDM;	
	MPI002B	BD connector (700077-001)				
	MPI002C	Abbott/Medix connector (700075-001)				
	MPI002D	Philips connector (22339301)				
	MPI002E	Edward/Baxter connector (2021197-001)				
	MPI002F	Medex connector				
	MPI002G	PVB connector				
	MPI002H	Medex Logical/MX960				

















Brand	Part No	Description	Picture	Compatibility	
Maquet	MPI027A	Utah connectr			
	MPI027B	BD connector			
	MPI027C	Abbott/Medix connector			
	MPI027D	Philips connector			
	MPI027E	Edward/Baxter connector			
	MPI027F	Medex connector			
	MPI027G	PVB connector			
	MPI027F	Medex Logical/MX960			
MEK	MPI026A	Utah connectr			
	MPI026B	BD connector (690-0021-00)			
	MPI026C	Abbott/Medix connector			
	MPI026D	Philips connector			
	MPI026E	Edward/Baxter connector			
	MPI026F	Medex connector			
	MPI026G	PVB connector			
	MPI026H	Medex Logical/MX960			














Brand	Part No	Description	Picture		Compatibility	
Mennen	MPI018A	Utah connectr			Mennen:Envoy, vitalogik 5500	
	MPI018B	BD connector				
	MPI018C	Abbott/Medix connector				
	MPI018D	Philips connector				
	MPI018E	Edward/Baxter connector				
	MPI018F	Medex connector				
	MPI018G	PVB connector				
	MPI018H	Medex Logical/MX960				
Mindray	MPI011A	Utah connectr			Mindray:PM6000,6800, T5, T8;Beneview T5/T6/T8 DPM 6, DPM 7	
	MPI011B	BD connector				
	MPI011C	Abbott/Medix connector				
	MPI011D	Philips connector				
	MPI011E	Edward/Baxter connector (0010-21-12179)				
	MPI011F	Medex connector				
	MPI011G	PVB connector				
	MPI011H	Medex Logical/MX960				






Brand	Part No	Description	Picture		Compatibility	
Nihon Kohden	MPI017A	Utah connectr			Nihon Kohden:Lifescopes 6/9/12/14 AP-311P -800PA BMS 2100-8500AOEC 3200-8120A/OMP-3221E -n7203/Mercury 9510K	
	MPI017B	BD connector (001C-30-70757)				
	MPI017C	Abbott/Medix connector				
	MPI017D	Philips connector				
	MPI017E	Edward/Baxter connector				
	MPI017F	Medex connector				
	MPI017G	PVB connector				
	MPI017H	Medex Logical/MX960				
Nihon Kohden	MPI015A	Utah connectr			Compatible with BSM2300 Series Lifescopes I/L, Other models need to be approved	
	MPI015B	BD connector				
	MPI015C	Abbott/Medix connector				
	MPI015D	Philips connector				
	MPI015E	Edward/Baxter connector				
	MPI015F	Medex connector				
	MPI015G	PVB connector				
	MPI015H	Medex Logical/MX960				

Brand	Part No	Description	Picture		Compatibility	
Philips /HP	MPI003A	Utah connectr			M1567A,M1568A;78352A, M1165A,Merlin Viridia 24,CMS 24,M3/M4,M1175A,M1350B,M1205A,M1350A,M1353A,8030A,8040A;HP:78205AD,78341,78342,78345,78346,78353,78354,783510,78551,78553,78801 78834,78833,Merlin 1300 series swithches automatically from 5uV to 40uV,M1006 to M1176A Press,M1006B,Hp,M P50/V24 Mindray: PM6000	
	MPI003B	BD connector				
	MPI003C	Abbott/Medix connector				
	MPI003D	Philips connector (M1634A)				
	MPI003E	Edward/Baxter connector				
	MPI003F	Medex connector				
	MPI003G	PVB connector				
	MPI003H	Medex Logical/MX960				
Siemens	MPI004A	Utah connectr			Siemens/Darger:SC6002XL,7000,SC6000,7000,9000, SC8000,SC9000XL	
	MPI004B	BD connector				
	MPI004C	Abbott/Medix connector				
	MPI004D	Philips connector				
	MPI004E	Edward/Baxter connector				
	MPI004F	Medex connector				
	MPI004G	PVB connector				
	MPI004H	Medex Logical/MX960				

Brand	Part No	Description	Picture		Compatibility	
Siemens	MPI005A	Utah connectr			Sirecust700, Sirecust630;Sirecust 400-1281, SC8000;Sirecust 400, series 961,1241,1280,1281	
	MPI005B	BD connector				
	MPI005C	Abbott/Medix connector				
	MPI005D	Philips connector				
	MPI005E	Edward/Baxter connector				
	MPI005F	Medex connector				
	MPI005G	PVB connector				
	MPI005H	Medex Logical/MX960				
Siemens	MPI008A	Utah connectr			Siemens	
	MPI008B	BD connector				
	MPI008C	Abbott/Medix connector				
	MPI008D	Philips connector				
	MPI008E	Edward/Baxter connector				
	MPI008F	Medex connector				
	MPI008G	PVB connector				
	MPI008H	Medex Logical/MX960				

Brand	Part No	Description	Picture	Compatibility
Spacelabs	MPI006A	Utah connectr		(Datex-Ohmeda TC-VTK.) For use with 690-0008-00, 690-0009-00,690-0010-00, 690-0011-00,690-0014-01, 690-0015-01 transducers and all Spacelabs. Use on 90496 requires adapter cable 700-0028-00.Mindray: PM7000,8000,9000,6201.; BCI,IVY,Ctitikon,Invivo,Me nnen,NDE,PB,Nellcor,Tekt ronix,Viatak;IVY:700,400 series,706,omni-trak4000,Ivnnvivo:Protocol system propeq 102,104,106,Acuity;Sieme ns:610,620;Tektronix,Vitat ek;Spacelabs:300/400/500 /600,Alpha7,Alpha,PC(903 00,90400,90600 series)PC express,700 surgical series,Tel/Vitatak,Scout TPC,Goldway:UT4000,Digi care All,Omni-Trak 4000;PHYSIO-CONTROL:VSM-1,VSM-4 VSM-5
	MPI006B	BD connector (690-0021-00)		
	MPI006C	Abbott/Medix connector		
	MPI006D	Philips connector		
	MPI006E	Edward/Baxter connector		
	MPI006F	Medex connector		
	MPI006G	PVB connector		
	MPI006H	Medex Logical/MX960		
Spacelabs	MPI014A	Utah connectr		Compatible: SpacelAb Ultraview 90496 , 91496 .
	MPI014B	BD connector		
	MPI014C	Abbott/Medix connector		
	MPI014D	Philips connector		
	MPI014E	Edward/Baxter connector		
	MPI014F	Medex connector		
	MPI014G	PVB connector		
	MPI014H	Medex Logical/MX960		

Brand	Part No	Description	Picture	Compatibility	
Marquette	MPI031A	GE-Marquette 2005772-001 IBP Convert cable		GE 11Pin-dual 11pin	
Marquette	MPI031B	GE 700-0028-00 IBP Convert cable		GE 11Pin-dual 6Pin	
Siemens	MPI032A	Siemens Draeger 5731281 IBP Convert Cable(13646)		Rec 16Pin-dual 10Pin	
Siemens	MP1032B	Siemens Draeger IBP Convert Cable		Rec 16Pin-dual 7Pin	
Siemens	MPI032C	Siemens IBP Convert Cable		10Pin-7Pin	
Siemens	MP1032D	Siemens Draeger 3368383E530U IBP Convert Cable		7Pin-10Pin	
Siemens	MPI032E	Siemens Draeger IBP Convert Cable		7Pin-Dual 10Pin	
Spacelabs	MPI033A	Spacelabs 700- 0028-00 IBP Convert cable		10Pin-Dual 6Pin	
Spacelabs	MPI033B	Spacelabs IBP Convert cable		10Pin-6Pin	
Mennen	MPI034A	Mennen 800-030- 270 IBP Convert cable		Round 10Pin-Dual 6Pin	
Mennen	MPI034B	Mennen IBP Convert cable		Round 10Pin-6Pin	
IBP Adaptor	MPI035	PVB-Abort			
IBP Adaptor	MPI036	PVB-Utah			

Brand	Part No	Description	Picture	Compatibility	
IBP Adaptor	MPI037	BD-Abort			
IBP Adaptor	MPI038	BD-Utah			
IBP Adaptor	MPI039	Phillips-Abort			
IBP Adaptor	MPI040	Phillips-Utah			
IBP Adaptor	MPI041	Utah-Abort			
IBP Adaptor	MPI042	Utah-Phillips	