



Medplus-Life plus,your best choice

Catalogue

Shenzhen Medplus Accessory Co.,Limited

Contact:Tina Wu,E-mail:info@spo2sale.com.Tel:+86-755-29093174


Address:2nd Floor,No.37,Ming'An Street,GongMing Town,

GuangMing New District,Shenzhen Of China


Catalogue

Spo2 Sensor

Nellcor


Nellcor	DS-100A MPS041A Nellcor adult finger clip spo2 sensor,3ft	
	MPS041B Nellcor adult soft tip spo2 sensor,3ft	
	D-YSPD MPS041C Nellcor pediatric finger clip spo2 sensor,3ft	
	MPS041D Nellcor pediatric soft tip spo2 sensor,3ft	
	D-YS MPS041E Nellcor neonate wrap spo2 sensor,3ft	
	D-YSE MPS041F Nellcor adult ear clip Spo2 sensor,3ft	
	MPS041G Nellcor animal tongue Spo2 sensor,3ft	

Compatibilities :Compatibility: All types, N40, N100, N290/5, N390, N3000, N6000

Nellcor	MPS042A Nellcor adult finger clip spo2 sensor,10ft	
	MPS042B Nellcor adult soft tip spo2 sensor,10ft	
	MPS042C Nellcor pediatric finger clip spo2 sensor,10ft	
	MPS042D Nellcor pediatric soft tip spo2 sensor,10ft	
	MPS042E Nellcor neonate wrap spo2 sensor,10ft	
	MPS042F Nellcor adult ear clip Spo2 sensor,10ft	
	MPS042G Nellcor animal tongue Spo2 sensor,10ft	

Compatibilities :Datascope:Accutorr Plus,Nellcor:NPB-290/5,N-390,N-395,N-3000,N6000,NPB-290-5,Colin BP 88/306

NellcorOximax

Nellcor Oximax	MPS043A Nellcor Oximax adult finger clip spo2 sensor,3ft	
	MPS043B Nellcor Oximax adult soft tip spo2 sensor,3ft	
	MPS043C Nellcor Oximax pediatric finger clip spo2 sensor,3ft	
	MPS043D Nellcor Oximax pediatric soft tip spo2 sensor,3ft	
	MPS043E Nellcor Oximax neonate wrap spo2 sensor,3ft	
	MPS043F Nellcor Oximax adult ear clip Spo2 sensor,3ft	
	MPS043G Nellcor Oximax animal tongue Spo2 sensor,3ft	

Compatibilities :GE-Medical System:Dinasmap,PRO1000;Nellcor:N550,N595,Require adapter EX023

Spo2 Adapter cable

	EC-8	MPE017A	<p>Artema: Diascope 2, Gamme 8000, Biosys BPM200, Ge-Marquette: dash, Eagle 3000, 4000, solar, TRAM X50-Modules, Hellige Eagle 1000, 2000, 3000 HP 78352A/C 78354A/C; Merlin: 1020A, 1024A, 1025A/B, 1205A, 1722A/B, 1732A/B, 2475B; Nellcor N100, N200, N180, N6000, N190, N20, N3000</p>		
Nellcor	DEC-8	MPE017B	<p>Datascope: Expert Datascope: Passport (Nellcor® SpO2); GE : Corometrics 188 Corometrics 129 Critikon Dinamap 9700 Critikon Dinamap 9710 Critikon Dinamap 9720 Plus Critikon Dinamap Compact Critikon Dinamap MPS Critikon Dinamap Pro 300 Critikon Dinamap Pro 400 Critikon Dinamap Select; IVY 405T, 405P; M.D.E.: Escort II OPT30 Prism; Nellcor: N100 N200 N180 N-20 NPB-40 NPB-75; Nonin: 8500; 8604; Welch Allyn Protocol: Propaq 206-EL, Protocol: Propaq 204-EL, Propaq 202-EL, Propaq 106EL</p>		
	DOC-10	MPE017C	<p>GE-MEDICAL SYSTEM: DINASMAP, Pro1000; NELLCOR : N550, N595 . Use with Nellcor Oximax tech sensor</p>		
Nellcor	SCP-10	MPE017D			